

IALUTEC® RED is a dietary supplement of pure hyaluronic acid with high molecular weight and red orange extract.

The pure hyaluronic acid, contained in IALUTEC® RED, is obtained by bacterial fermentation and not from animal sources. Hyaluronic acid is one of the fundamental components of human connective tissues; it is present in high concentrations in the synovial fluid, cartilage, skin and in the eyes. Experimental evidence shows that pure hyaluronic acid with high molecular weight, orally administered, reaches the peripheral tissues, the joints, the surface layers of the skin and especially the deep layers of the entire epidermis, where creams may find it harder to penetrate.¹ From 25 years of age onwards, human beings experience a progressive reduction in the production of hyaluronic acid due to physiological processes associated with aging. The red orange extract is a standardized extract obtained from the rind and the juice of three pigmented varieties of blood orange, Citrus Sinensis cultivated exclusively in a particular area surrounding the Etna volcano in Sicily (Italy). These oranges grow in an area particularly prone to extreme temperature and retain inside them a high antioxidant power as a defense from the external environment. The red orange extract contains a specific mixture titrated to 20% in active ingredients: polyphenols (anthocyanins, flavonoids, phenols) and ascorbic acid (or vitamin C). Ascorbic acid, also known as vitamin C, is an essential dietary element at the bottom of numerous cellular and tissue processes. In addition to bone and joint development and function, vitamin C is also the major factor involved in the synthesis of collagen, which together with hyaluronic acid, defines the structure for maintaining the elasticity of the connective tissue of blood vessels, the skin and muscles. Vitamin C is also recognized as a potent antioxidant and anti-UV.

DIRECTIONS For oral use only. Add 3 ml (36 drops) of IALUTEC® RED to a glass of water once a day for an initial period of 2 months. Thereafter add 2 ml (24 drops) once a day for the following months as maintenance. **GLUTEN FREE**

WARNINGS Keep out of reach of children. If you are pregnant, nursing or taking medications, consult your doctor before use. Do not use if seal is broken. Store at room temperature. Protect from moisture. Do not exceed recommended dose.

Net volume 50 ml / 1.69 fl oz

Brand product distributed in Italy by Eye Pharma SpA, Genoa

Distributed in the US by Eye Pharma Inc. 4770 Biscayne Blvd. Suite 400 Miami FL 33137. Made in Italy.

For more information visit www.eyepharma.com - info@eyepharma.com

1. Oe M, Mitsugi K, Odanaka W, Yoshida H, Matsuoka R, Seino S, Kanemitsu T, Masuda Y. Dietary hyaluronic acid migrates into the skin of rats. Scientific World Journal. 2014.